




DDX39A Antibody

Product Code	CSB-PA006620ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00148
Immunogen	Recombinant Human ATP-dependent RNA helicase DDX39A protein (178-427AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase DDX39A (EC 3.6.4.13) (DEAD box protein 39) (Nuclear RNA helicase URH49), DDX39A, DDX39
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DDX39A

Image



Western blot

All lanes: ATP-dependent RNA helicase

DDX39A antibody at 2μg/ml

Lane 1: A431 whole cell lysate

Lane 2: 293T whole cell lysate


Lane 3: HepG2 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 50, 37, 31 kDa

Observed band size: 50 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA006620ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA006620ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.