





EWSR1 Antibody

Product Code	CSB-PA007875ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01844
Immunogen	Recombinant Human RNA-binding protein EWS protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein EWS (EWS oncogene) (Ewing sarcoma breakpoint region 1 protein), EWSR1, EWS
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	EWSR1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA007875ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA007875ESR2HU at dilution of 1:100