





# FGF13 Antibody

<b>Product Code</b>	CSB-PA008619ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q92913
<b>Immunogen</b>	Recombinant Human Fibroblast growth factor 13 protein (1-245AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Rat
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Fibroblast growth factor 13 (FGF-13) (Fibroblast growth factor homologous factor 2) (FHF-2), FGF13, FHF2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	FGF13

## Image



Western blot  
All lanes: Fibroblast growth factor 13 antibody at 3µg/ml + Rat brain tissue  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 28, 22, 29, 23, 26 kDa  
Observed band size: 28 kDa



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA008619ESR1HU at dilution of 1:100

