





HDAC11 Antibody

Product Code	CSB-PA010237ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DB2
Immunogen	Recombinant Human Histone deacetylase 11 protein (1-347AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone deacetylase 11 (HD11) (EC 3.5.1.98), HDAC11
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HDAC11

Image



Western blot
 All lanes: HDAC11 antibody at 2µg/ml + Mouse brain tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 40, 34 kDa
 Observed band size: 40 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA010237ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.