



HS1BP3 Antibody

Product Code	CSB-PA703623ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53T59
Immunogen	Recombinant Human HCLS1-binding protein 3 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	HCLS1-binding protein 3 (HS1-binding protein 3) (HSP1BP-3), HS1BP3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HS1BP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.