





KMO Antibody

Product Code	CSB-PA012475ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15229
Immunogen	Recombinant Human Kynurenine 3-monooxygenase protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kynurenine 3-monooxygenase (EC 1.14.13.9) (Kynurenine 3-hydroxylase), KMO
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KMO

Image



Western blot
 All lanes: KMO antibody at 10µg/ml + Mouse heart tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 56, 55, 52 kDa
 Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA012475ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.