




LRG1 Antibody

Product Code	CSB-PA013082ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02750
Immunogen	Recombinant Human Leucine-rich alpha-2-glycoprotein protein (36-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich alpha-2-glycoprotein (LRG), LRG1, LRG
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LRG1

Image



Western blot
 All lanes: Leucine-rich alpha-2-glycoprotein antibody at 4µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 38 kDa
 Observed band size: 38 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.