





TARBP2 Antibody

Product Code	CSB-PA623007ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15633
Immunogen	Recombinant Human RISC-loading complex subunit TARBP2 protein (1-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RISC-loading complex subunit TARBP2 (TAR RNA-binding protein 2) (Trans-activation-responsive RNA-binding protein), TARBP2, TRBP
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	TARBP2

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA623007ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA623007ESR1HU at dilution of 1:100