





TGM2 Antibody

Product Code	CSB-PA023462ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21980
Immunogen	Recombinant Human Protein-glutamine gamma-glutamyltransferase 2 protein (438-687AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-glutamine gamma-glutamyltransferase 2 (EC 2.3.2.13) (Tissue transglutaminase) (Transglutaminase C) (TG(C)) (TGC) (TGase C) (Transglutaminase H) (TGase H) (Transglutaminase-2) (TGase-2), TGM2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TGM2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA023462ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA023462ESR1HU at dilution of 1:100

