




XRCC4 Antibody

Product Code	CSB-PA614413ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13426
Immunogen	Recombinant Human DNA repair protein XRCC4 protein (1-336AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA repair protein XRCC4 (X-ray repair cross-complementing protein 4), XRCC4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	XRCC4

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA614413ESR2HU at dilution of 1:100