





ZEB2 Antibody

Product Code	CSB-PA026425DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60315
Immunogen	Recombinant Human Zinc finger E-box-binding homeobox 2 protein (1050-1214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger E-box-binding homeobox 2 (Smad-interacting protein 1) (SMADIP1) (Zinc finger homeobox protein 1b), ZEB2, KIAA0569 SIP1 ZFH1B ZFX1B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ZEB2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA026425DSR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA026425DSR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.