




ACOT7 Antibody

Product Code	CSB-PA001169ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00154
Immunogen	Recombinant Human Cytosolic acyl coenzyme A thioester hydrolase protein (151-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic acyl coenzyme A thioester hydrolase (EC 3.1.2.2) (Acyl-CoA thioesterase 7) (Brain acyl-CoA hydrolase) (BACH) (CTE-IIa) (CTE-II) (Long chain acyl-CoA thioester hydrolase), ACOT7, BACH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACOT7

Image


Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA001169ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.