





ACSM3 Antibody

Product Code	CSB-PA709601ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53FZ2
Immunogen	Recombinant Human Acyl-coenzyme A synthetase ACSM3, mitochondrial protein (387-586AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl-coenzyme A synthetase ACSM3, mitochondrial (EC 6.2.1.2) (Acyl-CoA synthetase medium-chain family member 3) (Butyrate--CoA ligase 3) (Butyryl-coenzyme A synthetase 3) (Middle-chain acyl-CoA synthetase 3) (Protein SA homolog), ACSM3, SAH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACSM3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA709601ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA709601ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.