



AMBRA1 Antibody

Product Code	CSB-PA883655ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0C7
Immunogen	Recombinant Human Activating molecule in BECN1-regulated autophagy protein 1 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Activating molecule in BECN1-regulated autophagy protein 1, AMBRA1, KIAA1736
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AMBRA1