





BCAS2 Antibody

Product Code	CSB-PA002597ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75934
Immunogen	Recombinant Human Pre-mRNA-splicing factor SPF27 protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-mRNA-splicing factor SPF27 (Breast carcinoma-amplified sequence 2) (DNA amplified in mammary carcinoma 1 protein) (Spliceosome-associated protein SPF 27), BCAS2, DAM1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	BCAS2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA002597ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA002597ESR2HU at dilution of 1:100