




CD200R1 Antibody


Product Code	CSB-PA851547ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD46
Immunogen	Recombinant Human Cell surface glycoprotein CD200 receptor 1 protein (226-325AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface glycoprotein CD200 receptor 1 (CD200 cell surface glycoprotein receptor) (Cell surface glycoprotein OX2 receptor 1), CD200R1, CD200R CRTR2 MOX2R OX2R
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD200R1

Image



Western blot

All lanes: CD200R1 antibody at 6µg/ml + Mouse lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 37, 22, 19, 40 kDa
 Observed band size: 37 kDa



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA851547ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.