





DNAJC2 Antibody

Product Code	CSB-PA858715ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99543
Immunogen	Recombinant Human DnaJ homolog subfamily C member 2 protein (1-140AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 2 (M-phase phosphoprotein 11) (Zuotin-related factor 1) [Cleaved into: DnaJ homolog subfamily C member 2, N-terminally processed], DNAJC2, MPHOSPH11 MPP11 ZRF1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DNAJC2

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA858715ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA858715ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.