



HPSE2 Antibody

Product Code	CSB-PA823898ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWQ2
Immunogen	Recombinant Human Inactive heparanase-2 protein (42-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive heparanase-2 (Hpa2), HPSE2, HPA2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HPSE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.