




PARL Antibody

Product Code	CSB-PA880960ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H300
Immunogen	Recombinant Human Presenilins-associated rhomboid-like protein, mitochondrial protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Presenilins-associated rhomboid-like protein, mitochondrial (EC 3.4.21.105) (Mitochondrial intramembrane cleaving protease PARL) [Cleaved into: P-beta (Pbeta)], PARL, PSARL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PARL

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA880960ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.