



RASGRP3 Antibody

Product Code	CSB-PA808539ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IV61
Immunogen	Recombinant Human Ras guanyl-releasing protein 3 protein (541-690AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras guanyl-releasing protein 3 (Calcium and DAG-regulated guanine nucleotide exchange factor III) (Guanine nucleotide exchange factor for Rap1), RASGRP3, GRP3 KIAA0846
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASGRP3