



Pectate lyase 1 Antibody

Product Code	CSB-PA856862LA01CHAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96385
Immunogen	Recombinant Chamaecyparis obtusa Pectate lyase 1 protein (22-354AA)
Raised In	Rabbit
Species Reactivity	Chamaecyparis obtusa
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pectate lyase 1 (EC 4.2.2.2) (Major pollen allergen Cha o 1) (allergen Cha o 1)
Immunogen Species	Chamaecyparis obtusa (Hinoki false-cypress)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.