



IFNAR2 Antibody, Biotin conjugated

Product Code	CSB-PA11439D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48551
Immunogen	Recombinant Human Interferon alpha/beta receptor 2 protein (27-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha/beta receptor 2 (IFN-R-2) (IFN-alpha binding protein) (IFN-alpha/beta receptor 2) (Interferon alpha binding protein) (Type I interferon receptor 2), IFNAR2, IFNABR IFNARB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNAR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.