



Batf3 Antibody

Product Code	CSB-PA002572LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97876
Immunogen	Recombinant Rat Basic leucine zipper transcriptional factor ATF-like 3 protein (1-133AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Basic leucine zipper transcriptional factor ATF-like 3 (B-ATF-3) (Jun dimerization protein 1) (JDP-1), Batf3, Jdp1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Epigenetics and Nuclear Signaling
Target Names	Batf3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.