



stxB2 Antibody

Product Code	CSB-PA357663HA01ECC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09386
Immunogen	Recombinant Enterobacteria phage 933W Shiga-like toxin 2 subunit B protein (20-89AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage 933W
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Shiga-like toxin 2 subunit B (SLT-2 B subunit) (SLT-2b) (SLT-IIb) (Verocytotoxin 2 subunit B) (Verotoxin 2 subunit B), stxB2, stx2B
Immunogen Species	Enterobacteria phage 933W (Bacteriophage 933W)
Research Area	Others
Target Names	stxB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.