



# Minpp1 Antibody

<b>Product Code</b>	CSB-PA013831LA01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O35217
<b>Immunogen</b>	Recombinant Rat Multiple inositol polyphosphate phosphatase 1 protein (31-481AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Multiple inositol polyphosphate phosphatase 1 (EC 3.1.3.62) (2,3-bisphosphoglycerate 3-phosphatase) (2,3-BPG phosphatase) (EC 3.1.3.80) (Inositol (1,3,4,5)-tetrakisphosphate 3-phosphatase) (Ins(1,3,4,5)P(4) 3-phosphatase), Minpp1, Mipp1
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	Minpp1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.