





LYZ Antibody

Product Code	CSB-PA02769A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61626
Immunogen	Recombinant Human Lysozyme C protein (19-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysozyme C (EC 3.2.1.17) (1,4-beta-N-acetylmuramidase C), LYZ, LZM
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LYZ

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA02769A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA02769A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.