





LRRN1 Antibody

Product Code	CSB-PA747638LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXK5
Immunogen	Recombinant Human Leucine-rich repeat neuronal protein 1 protein (26-631AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat neuronal protein 1 (Neuronal leucine-rich repeat protein 1) (NLRR-1), LRRN1, KIAA1497
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRRN1

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA747638LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA747638LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.