




# CEP63 Antibody

<b>Product Code</b>	CSB-PA839376LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96MT8
<b>Immunogen</b>	Recombinant Human Centrosomal protein of 63 kDa protein (1-495AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Centrosomal protein of 63 kDa (Cep63), CEP63
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	CEP63

## Image



### Western blot

All lanes: CEP63 antibody at 12µg/ml

Lane 1: K562 whole cell lysate

Lane 2: HeLa whole cell lysate

Lane 3: NIH/3T3 whole cell lysate


Lane 4: Raji whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 82, 64, 58, 56 kDa

Observed band size: 64 kDa



### Immunohistochemistry of paraffin-embedded

human prostate cancer using CSB-

PA839376LA01HU at dilution of 1:100



**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.