



Serpinf1 Antibody

Product Code	CSB-PA021084DA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q80ZA3
Immunogen	Recombinant Rat Alpha-2 antiplasmin protein (20-418AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2 antiplasmin, Protein Serpinf1, Serine (Or cysteine) peptidase inhibitor, clade F, member 1, Serine (Or cysteine) peptidase inhibitor, clade F, member 1, isoform CRA_c, Serpinf1, Dmrs91, rCG_34442
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Serpinf1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.