




# SLC27A6 Antibody

<b>Product Code</b>	CSB-PA897533ESR2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y2P4
<b>Immunogen</b>	Recombinant Human Long-chain fatty acid transport protein 6 protein (370-619AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Long-chain fatty acid transport protein 6 (FATP-6) (Fatty acid transport protein 6) (Fatty-acid-coenzyme A ligase, very long-chain 2) (Solute carrier family 27 member 6) (Very long-chain acyl-CoA synthetase homolog 1) (VLCSH1) (hVLCs-H1), SLC27A6, ACSVL2 FACVL2 FATP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SLC27A6

## Image



### Western blot

All lanes: SLC27A6 antibody at 4µg/ml + MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 70 kDa

Observed band size: 70 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.