




SLC27A6 Antibody

Product Code	CSB-PA897533ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2P4
Immunogen	Recombinant Human Long-chain fatty acid transport protein 6 protein (370-619AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain fatty acid transport protein 6 (FATP-6) (Fatty acid transport protein 6) (Fatty-acid-coenzyme A ligase, very long-chain 2) (Solute carrier family 27 member 6) (Very long-chain acyl-CoA synthetase homolog 1) (VLCSH1) (hVLCs-H1), SLC27A6, ACSVL2 FACVL2 FATP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC27A6

Image



Western blot

All lanes: SLC27A6 antibody at 4µg/ml + MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 70 kDa
 Observed band size: 70 kDa