




SPAG5 Antibody

Product Code	CSB-PA853491ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96R06
Immunogen	Recombinant Human Sperm-associated antigen 5 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm-associated antigen 5 (Astrin) (Deepest) (Mitotic spindle-associated protein p126) (MAP126), SPAG5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPAG5

Image



Western blot


All lanes: SPAG5 antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 134 kDa

Observed band size: 134 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA853491ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.