



# NDEL1 Antibody

<b>Product Code</b>	CSB-PA872414ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9GZM8
<b>Immunogen</b>	Recombinant Human Nuclear distribution protein nudE-like 1 protein (1-345AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Nuclear distribution protein nudE-like 1 (Protein Nudel) (Mitosis-associated protein 1), NDEL1, EOPA MITAP1 NUDEL
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	NDEL1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.