





PAX7 Antibody

Product Code	CSB-PA017493ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23759
Immunogen	Recombinant Human Paired box protein Pax-7 protein (291-520AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired box protein Pax-7 (HuP1), PAX7, HUP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PAX7

Image



Western blot
All lanes: PAX7 antibody at 4μg/ml + 293T whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 56, 57 kDa
Observed band size: 56 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA017493ESR1HU at dilution of 1:100