



GMPR2 Antibody

Product Code	CSB-PA873740LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2T1
Immunogen	Recombinant Human GMP reductase 2 protein (1-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GMP reductase 2 (GMPR 2) (EC 1.7.1.7) (Guanosine 5'-monophosphate oxidoreductase 2) (Guanosine monophosphate reductase 2), GMPR2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GMPR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.