



GMPR2 Antibody

Product Code	CSB-PA873740LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2T1
Immunogen	Recombinant Human GMP reductase 2 protein (1-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GMP reductase 2 (GMPR 2) (EC 1.7.1.7) (Guanosine 5'-monophosphate oxidoreductase 2) (Guanosine monophosphate reductase 2), GMPR2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GMPR2