





# CD38 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA004929ESR1HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                          |
| <b>Uniprot No.</b>         | P28907   |
| <b>Immunogen</b>           | Recombinant Human ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 protein (75-300AA) |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA, IHC; Recommended dilution: IHC:1:20-1:200                                       |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                      |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Immunology   |
| <b>Target Names</b>        | CD38   |

## Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA004929ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA004929ESR1HU at dilution of 1:100