



# The metabolite of furaltadone Monoclonal Antibody

<b>Product Code</b>	CSB-MA007111I0m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Tested Applications</b>	ELISA
<b>Form</b>	liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95% , Protein G purified
<b>Isotype</b>	IgG1
<b>Clonality</b>	monoclonal
<b>Alias</b>	CPAMOZ, The metabolite of furaltadone
<b>Product Type</b>	Small Molecules Antibodies
<b>Target Names</b>	AMOZ
<b>Accession NO.</b>	1C1G8