



The metabolite of furaltadone Monoclonal Antibody

Product Code	CSB-MA007111I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 0.01M PBS, PH 7.4
Purification Method	>95% , Protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	CPAM0Z, The metabolite of furaltadone
Product Type	Small Molecules Antibodies
Target Names	AM0Z
Accession NO.	1C1G8