



Tmprss4 Antibody

Product Code	CSB-PA844988LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VCA5
Immunogen	Recombinant Mouse Transmembrane protease serine 4 protein (52-435AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protease serine 4 (EC 3.4.21.-) (Channel-activating protease 2) (mCAP2), Tmprss4, Cap2
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Tmprss4