




SYT11 Antibody

Product Code	CSB-PA883595ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BT88
Immunogen	Recombinant Human Synaptotagmin-11 protein (37-200AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptotagmin-11 (Synaptotagmin XI) (SytXI), SYT11, KIAA0080
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SYT11

Image



Western blot

All lanes: SYT11 antibody at 14µg/ml

Lane 1: COLO205 whole cell lysate

Lane 2: Mouse skeletal muscle tissue

Lane 3: MCF-7 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.