




TNNC2 Antibody

Product Code	CSB-PA024010ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02585
Immunogen	Recombinant Human Troponin C, skeletal muscle protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Troponin C, skeletal muscle, TNNC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNNC2

Image



Western blot

All lanes: TNNC2 antibody at 10µg/ml + Mouse skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 18 kDa

Observed band size: 18 kDa