





TNNC2 Antibody

Product Code	CSB-PA024010ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02585
Immunogen	Recombinant Human Troponin C, skeletal muscle protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Troponin C, skeletal muscle, TNNC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNNC2

Image



Western blot
 All lanes: TNNC2 antibody at 10µg/ml + Mouse skeletal muscle tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 18 kDa
 Observed band size: 18 kDa



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA024010ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.