




SLC25A15 Antibody

Product Code	CSB-PA897578ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y619
Immunogen	Recombinant Human Mitochondrial ornithine transporter 1 protein (1-301AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial ornithine transporter 1 (Solute carrier family 25 member 15), SLC25A15, ORNT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC25A15

Image



Western blot

All lanes: SLC25A15 antibody at 6μg/ml

Lane 1: Mouse liver tissue


Lane 2: Mouse stomach tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa



Immunohistochemistry of paraffin-embedded

human placenta tissue using CSB-PA897578ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.