



Recombinant Human Proliferation marker protein Ki-67 (MKI67), partial

Product Code	CSB-EP014597HU(A5)
Relevance	Thought to be required for maintaining cell proliferation.
Abbreviation	Recombinant Human MKI67 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46013
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	TKRHLRTRVQKVQVKEEPSAVKFTQTSGETTDADKEPAGEDKGIKALKESAKQ TPAPAASVTGSRRRPRAPRESAQAIEDLAGFKDPAAGHTEESMTDDKTTKIPC KSSPELEDTATSSKRRPRTRAQKVEVKEELLAVGKLTQTSGETTHTDKEPVGE GKGTKAFKQPAKRKLD AEDVIGSRRQPRAPKEKAQPLEDLASFQELSQTTPGH TEELANGAADSFTSAPKQTPDSGKPLKISRRVLRAPKVEPVGDVVSTRDPVKS QSKSNTSLPPLPFKRGGGKDGSVTGKRLRCMPAPEEIVEELPASKKQRVAP RARGKSSEPVVIMKRSLRTSAKRIEPAEELNSNDMKTNKEEHKLQDSVPENKG ISLRSRRQNKTEAEQQITEVFVLAERIEINRNEKKPMKTSPEMDIQNPDDGARK PIPRDKVTENKRCLRSARQNESSQPKVAEESGGQKSAKVLMMQNQKKGKGEAG NSDSMCLRSRKTQSQAASTLESKSVQRVTRSVKRC AENPKKAEDNVCVKKI RTRSHRDSEDI
Research Area	Cancer
Source	E.coli
Target Names	MKI67
Expression Region	2720-3256aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	72.1 kDa
Protein Length	Partial
Image	

