



AGA2 Antibody

Product Code	CSB-PA335766HA01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32781
Immunogen	Recombinant Saccharomyces cerevisiae A-agglutinin-binding subunit protein (19-87AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	A-agglutinin-binding subunit, AGA2
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Others
Target Names	AGA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.