





GTF2H2C Antibody

Product Code	CSB-PA757553ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1K8
Immunogen	Recombinant Human General transcription factor IIH subunit 2-like protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	General transcription factor IIH subunit 2-like protein (General transcription factor IIH polypeptide 2-like protein), GTF2H2C; GTF2H2C_2GTF2H2D
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GTF2H2C

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA757553ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA757553ESR2HU at dilution of 1:100