




STK38 Antibody

Product Code	CSB-PA614885ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15208
Immunogen	Recombinant Human Serine/threonine-protein kinase 38 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 38 (EC 2.7.11.1) (NDR1 protein kinase) (Nuclear Dbf2-related kinase 1), STK38, NDR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STK38

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA614885ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.