



# UBE2D4 Antibody

<b>Product Code</b>	CSB-PA896498DSR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y2X8
<b>Immunogen</b>	Recombinant Human Ubiquitin-conjugating enzyme E2 D4 protein (1-147AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ubiquitin-conjugating enzyme E2 D4 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme D4) (HBUCE1) (Ubiquitin carrier protein D4) (Ubiquitin-protein ligase D4), UBE2D4, UBCH5D
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	UBE2D4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.