



SLC27A2 Antibody

Product Code	CSB-PA021534ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14975
Immunogen	Recombinant Human Very long-chain acyl-CoA synthetase protein (30-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Very long-chain acyl-CoA synthetase (VLACS) (VLCS) (EC 6.2.1.-) (Fatty acid transport protein 2) (FATP-2) (Fatty-acid-coenzyme A ligase, very long-chain 1) (Long-chain-fatty-acid--CoA ligase) (EC 6.2.1.3) (Solute carrier family 27 member 2) (THCA-CoA ligase) (Very long-chain-fatty-acid-CoA ligase), SLC27A2, ACSVL1 FACVL1 FATP2 VLACS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC27A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.